




THD LIGHTSCOPE RECTOSCOPE LED



Use

The rectoscope THD LightScope LED is a DISPOSABLE medical device, indicated for analysis of ano-rectal canal and the distal region of the sigmoid.

General Description

The rectoscope THD LightScope LED combines into a single disposable device a hand lens and an integrated

powerful LED light source, in order to have a clear and detailed vision of inspection of the area.

The THD Light-Scope rectoscope is ready for use and is also available with the blowing pump included in the sales package.

Features

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- Self-light, with an integrated light source
- 4x magnifier lens
- Two length available 200 mm and 250 mm (7.87 in and 9.84 in)
- Available also with disposable insufflator
- Ready to use
- Light (55 g/0.12 pound)
- FDA Regulatory Class II
- Individually packaged
- Latex Free



Technical Data

P/N. 800097
 P/N. 800098
 P/N. 800099
 P/N. 800100

US

Advantages

- It doesn't require an additional light source
- It doesn't require a preparation before its use
- It delete the risk of cross-contamination

Dimensional Features		
THD LIGHTSCOPE RECTO LED		
	Mini	Maxi
Diameter (B)	23 mm (0.90 in)	23 mm (0.90 in)
Length (A)	200 mm (7.87 in)	250 mm (9.84 in)

THD LightScope Recto LED Sterile		
Code	Qty	Denomination
800097 - 10	10	THD LightScope Recto Maxi LED
800098 - 10	10	THD LightScope Recto Mini LED
800099 - 10	10	THD LightScope Recto Maxi LED & bulb
800100 - 10	10	THD LightScope Recto Mini LED & bulb