




## THD LIGHTSCOPE RECTOSCOPE LED



### Use

The rectoscope THD LightScope LED is a DISPOSABLE medical device, indicated for analysis of ano-rectal canal and the distal region of the sigmoid.

### General Description

The rectoscope THD LightScope LED combines into a single disposable device a hand lens and an integrated

powerful LED light source, in order to have a clear and detailed vision of inspection of the area.

The THD Light-Scope rectoscope is ready for use and is also available with the blowing pump included in the sales package.

### Features

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- Self-light, with an integrated light source
- 4x magnifier lens
- Two length available 200 mm and 250 mm (7.87 in and 9.84 in)
- Available also with disposable insufflator
- Ready to use
- Light (55 g/0.12 pound)
- FDA Regulatory Class II
- Individually packaged
- Latex Free



# Technical Data

P/N. 800097  
 P/N. 800098  
 P/N. 800099  
 P/N. 800100

US

## Advantages

- It doesn't require an additional light source
- It doesn't require a preparation before its use
- It delete the risk of cross-contamination

Dimensional Features		
THD LIGHTSCOPE RECTO LED		
	<b>Mini</b>	<b>Maxi</b>
Diameter (B)	23 mm (0.90 in)	23 mm (0.90 in)
Length (A)	200 mm (7.87 in)	250 mm (9.84 in)

THD LightScope Recto LED		
Code	Qty	Trade Name
800097	25	THD LS MAXI-RECTO (LED NS)
800098	10	THD LS MINI-RECTO (LED NS)
800099	10	THD LS MAXI-RECTO (LED INS/NS)
	25	
	50	
800100	25	THD LS MINI-RECTO (LED INS/NS)