




THD LIGHTSCOPE LED ANOSCOPE



Use

The LightScope LED anoscope is a DISPOSABLE medical device deemed for investigation of the anal canal.

General Description

The LightScope LED anoscope combines in a single device a light anoscope and a powerful, fully integrated, LED light source, in order to have a clear and detailed vision of the inspection field.

The LightScope LED anoscope is ready to use and does not need for other tools
Available with truncated or beveled tip.

Features

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- **Self-light, with a permanent integrated light source**
- Two ends available (truncated and bevelled)
- Ready to use
- Light (40 g/0.09 pound)
- FDA Regulatory Class II
- Individually packaged
- Latex Free



Advantages

- It is not need for an external light source
- No need for set-up before use
- It avoids the risk of cross-contamination

Dimensional Features	
THD LIGHTSCOPE LED ANOSCOPE	
Operative length (A)	70 mm (2.76 in)
Outer diameter (B)	25 mm (0.98 in)
Internal diameter (C)	18 mm (0.71 in)
total length (D)	110 mm (4.33 in)
Bevelled total length (Af)	55 mm (2.17 in)

THD LightScope LED Anoscope			
Code	Qty	Trade Name	Denomination
800092 - 10	10	THD LIGHT-SCOPE ANO (LED TR/NS)	THD LightScope Ano Truncated LED
800092 - 25	25		
800092 - 50	50		
800093 - 10	10	THD LIGHT-SCOPE ANO (LED BE/NS)	THD LightScope Ano Bevelled LED
800093 - 25	25		
800093 - 50	50		