



# THD LIGHTSCOPE LED ANOSCOPE STERILE



## Use

The LightScope LED anoscope is a DISPOSABLE and STERILE medical device deemed for investigation of the anal canal.

## General Description

The LightScope LED anoscope combines in a single device a light anoscope and a powerful, fully integrated, LED light source, in order to have a clear and detailed vision of the inspection field.

The LightScope LED anoscope is ready to use and does not need for other tools  
Available with truncated or beveled tip.

## Features

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- Self-light, with an integrated light source
- Two ends available (truncated and beveled)
- Ready to use
- Light (40 g/0.09 pound)
- FDA Regulatory Class II
- Individually packaged
- Latex Free
- Sterilized by Ethylene Oxide



Advantages

- It is not need for an external light source
- No need for set-up before use
- It avoids the risk of cross-contamination
- Sterile Device

Dimensional Features	
THD LIGHTSCOPE LED ANOSCOPE STERILE	
Diameter (B)	25 mm (0.98 in)
Length (A)	70 mm (2.76 in)

THD LightScope LED Anoscope Sterile		
TRUNCATED		BEVELLED
Code	Qty	Denomination
800112 - 10	10	THD LightScope Ano Truncated LED Sterile
800113 - 10	10	THD LightScope Ano Bevelled LED Sterile

THD360502\_US\_Rev00\_20190424\_THD\_Lightscope\_Ano\_LED\_Sterile.docx