




THD SUTURE PACK (6 sutures)



Use

The THD Suture is used for the treatment of the haemorrhoidal disease, using the Transanal Haemorrhoidal Dearterialisation technique.

General Description

THD Sutures with needle is are absorbable surgical sutures, U.S.P.

Braided and coated Polyglycolic acid.

Dimensional Features	
Wire Dimension	(2-0) U.S.P. 3 Metric
Needle	5/8 Circle Round Bodied 26.5 mm
Total Length	700 mm (≈ 27 in)

THD SUTURE			
Qty	P/N	Description	510(k)
1	880000	THD Suture pack (6 sutures)	K041515