



ROLLING CART FOR THD REVOLUTION



Use

The Rolling Cart creates a mobile stand during procedure and storage of the THD Revolution (doppler dearterialiser unit).

General Description

The Rolling Cart creates greater mobility, ease of orientation and negates the need to move other instruments in the room.

The cart is equipped with a handle for better control and ease of use.

Dimensional Features	
Dimension	Height: 1000 mm (39 in); Pivoting diameter 610 mm (24 In)
Weight	13 kg (28.7 lb)

THD SUTURE		
Qty	P/N	Description
1	TUSX88014	Rolling Cart for THD Revolution